Leica STELLARIS 8 Confocal/FLIM/tauSTED Operating Instructions

System Overview (Confocal Application)



Hardware information





1	White Light Laser (WLL)		External light source		
2	Scan Head	12	Workstation		
3	Microscope	14	External cooling (chiller)		
4	Monitor(s)	15	Supply Unit		
5	Control panel				

Sample mounting and fluorescence observation through eyepieces





- 1. Tilt the illumination column.
- 2. Verify that the correct sample holder stage inset is correctly secured in the XY stage.





3. If using immersion lens, clean the lens and apply proper immersion.

4. Never clean dry lens!

5. If suspected dirty, or if immersion gets onto a dry lens, contact KCCI staff immediately!

Sample mounting and fluorescence observation through eyepieces

6. Secure your sample, thin coverslip Facing Down. The objectives cannot focus through a glass slide. If using immersion lens, quickly bring the immersion to contact with the sample using coarse/fast Z-focus (and XY-stage) adjustment on the "Smart Move" controller.



Coarse/fine XY-stage

7. Using the Aperture (aperture diaphragm) or Field (field diaphragm) keys, the motorized diaphragms can be changed at any time. The display on the Touch Screen changes accordingly.



Coarse/fine /Z-Focus



System Start-up



- Turn on the power (wait for few seconds) and laser button as shown in the figure 1.
 *Note: Do not touch the key.
- 2. Turn ON the computer and wait.
 3. Log into the computer with your own KCCI user account.
 4. Record the starting time at the KCCI website



5. Start the system by double clicking the LAS X icon on the desktop.

*Note: The LAS X software is used to control all system functions and acts as the link to the individual hardware components. Image acquisition, image analysis and image processing are carried out using LAS X.

	Peira
Leica Application Suite X 4.5.0.25531	MICROSYSTEMS
Configuration : Microscope :	machine.xlhw 🗘 DMI8 🗘
STED : Load settings at startup :	COFF.
F	ОК Cancel

Leica Application Suite X 4.5.0.25531 Initializing SuperContVisible Light... Initializing ULL... Initializing WLL... Initializing COM11... Initializing DrvV24Device 0...

6. Select the Machine as configuration and DM18 in the Microscope section as shown.

7. Always OFF for the confocal usage.

8. Now start LAS X by clicking the OK button.

9. Wait Until the software initialization completes .



10. Select No, if you are NOT using for tiling, etc. For Tiling you should initialize the stage.

11. Wait for the LAS X Graphical User Interface window.

Using the Microscope Touchscreen



On touch screen click the highlighted icon for basic microscope settings. Select the status/illumination menu using light source icon. Adjust the field diaphragms using + and – keys. FIM: Fluorescence intensity manager.

FINI: Fluorescence intensity manager.

IL: Incident light (Fluo: Fluorescence light)

TL: Transmitted light

Contrast Tab or Modality Tab:

Transmitted light: Bright field (BF) **Incident light**: Fluorescence (FLUO) **Combination**: Fluorescence-Phase (FLUO-PH) and (CS)

The shutter can be operated using the buttons on the touch screen.

Click the FLUO for fluorescence incident light method

The available filter cubes are displayed. Select the desired cube using the corresponding button.



To see on the screen, use these settings. Combi Incident Transmitted FLUO-PH FLUO CS Gray: Shutter closed Yellow: Shutter open xZ У IL-Shutter **TL-Shutter** FLUO-Filtercubes EMP_BF LED 405 GFP RHOD LP Leica





Using the microscope touchscreen



Objective Tab

- 1. Shows the total magnification on the eyepieces.
- 2. Available objectives: Here 5.

Click the objective button from the touch screen and select the required objective.

- Use only <u>Type F immersion liquid</u> for the objective if required.
- 3. Next part of the panel shows whether the light is going to the eyepieces or to the camera.





Structure of the LAS X Graphical User Interface



Click the <u>configuration button</u> to activate the configuration



Click the Hardware button to select the Panning, Resolution etc.

Jein -0-	STELLARIS 8 🗘 Configuration Acquire	Process Quantify Analysis 🖯
IPS Beam Path	🔅 Hardware Setti	ngs
Laser Config	▼ Panning Step Size	 Data Transfer Mode Direct Direct Overflow Maximum Integration Time
Objective Super-Z Hardware	✓ Resolution Bit Depth ✓ Z-Movement Activate Z-Movement Limitation	 Online Maximum Projection during Acquisition Online Maximum Projection during Acquisition Online Maximum Projection during Acquisition Image Orientation Flip X
User Config User Config Dye Database	 Bidirectional XYZT Acquisition Manual Microscope Control Enable during Acquisition of Series Detector Auto Selection 	 LDM/Wizard Stop Behaviour Wait until frame is ready

Click the **Objective Icon** to select the desired objective the configuration panel.



To select the required fluorophore, click the Dye Database.



Click the Laser configuration icon and turn on the white light laser (WLL) from the dialogue box appeared on the screen.

<u>Peter</u>	-0	7/	STELLARIS 8 🗘	Configuration	Acquire	Process	Quantify	Analysis 🖨
IPS		Beam Path						
	*							
Laser C	onfig	USB Panel						
Objec	tive	Super-Z						
Q	×							
Hardv	vare							
			🔹 Gurrenti	y available	e Lasers			
Ť			▼ Adjust Laser Settings	Ì				*
User Co	onfig	Memory						
		••	405 Diode:	DFF				
Dve Dat	abase	CAMServer	WLL: ON				o	85.00 %
Uye Dat	dease	GAMBERET						

Click the Acquire panel to image the specimen. The Acquire panel is shown in the figure.



Configuration	Acquire	Process	Quantify	, >	Analysis 🗎
		Z-stack imaging		Time	-lapse
Open projectsAcquisitionV Acquis tion ModeImage SizexyztFLMImage SizeV XY: 512x512 400 Hz 1.00 1.00 AUImage SizeV XY: 512x512 400 Hz 1.00 1.00 AUImage SizeSpeed:400 + Image SizeY XY: 512x512 400 Hz 1.00 1.00 AUImage SizeSpeed:400 + Image SizeY XY: 512x512 400 Hz 1.00 1.00 AUImage SizeY Y: 512x512 5 µm x 581.25 µm x 581.25 µm x 581.25 µmY Y: 512x512 5 µm x 581.25 µm x 58	 Choose the acquisition mode from the dropdown list. Use the XY panel to set the image format, zoom, line averaging, pinhole. 	 ▼ Z-Stack : Begin ← End ← C Position C Position C Position C Z Size [O.0 Z Size [O.0 Re-Cer 	① ↓ ① ↓ ① ↓ ① ↓ ① ↓ ① ↓ ① ↓ ① ① ↓ ① ① ① ① ① ① ① ① ①	 ▼ t: 1 00:00:01.293 Time Interval Minimize Acquire Until Stop Duration ≈ Stacks 	h 00:00:01.293 h 0: 0: 0: 1:293 DAYS HOURS MIN SEC 1/1000 0: 0: 0: 1:293 DAYS HOURS MIN SEC 1/1000 1 ◆

Loading your imaging method



• Click the Acquisition panel under Acquire and click the Live to check the specimen or Capture Image button to record the image.

- NOTE: The method must be set up by the KCCI staff and you must be trained by the KCCI staff to use the method.
- The easiest way of loading a method is to open a previously saved image ('lif') under the Acquire 'Open Projects" and then click the 'Apply' button to load the imaging method.





Acquire panel: Selection of the detectors, Dye, and Excitation wavelength.



- To select the fluorophore from the available top ribbon, hold the required dye button and drag it to Setting 1 (bottom panel).
- The particular hybrid detector will appear at the bottom of the selected fluorophore.

Acquire panel: Selection of the detectors, Dye, and Excitation wavelength.



- To remove the unwanted fluorophore active the detector ribbon (below the fluorophore) from the setting 1 panel and click the delete button.
- To remove the unwanted excitation wavelength, click the excitation line from the setting 1 panel and click the delete button.

Multi-fluorophores: Sequential Imaging Method



Bright Field



• Turn ON the Trans PMT.

*Note: Do not select the dye in the case of Bright Field.

- Drag the Add Laser button to the setting 1 panel to select the excitation light. Make sure the white light laser is ON.
- Click Live to view the specimen and click Capture to record the image.

Take a Z-Stack

*Note: Z-stack is different from Z-galvo

Acquisition Open projects 1 * Acquisition Mode xvz 🗘 Stack : 40µm | 21 Steps 3 Begin : [µm] 24.79 End : [µm] -15.20 Z Position [µm] -15.10 Beain-Z Size [µm] 40.00 Re-Center Stack Direction (Z) z-galvo 🖨 + **0** 🕀 Number of Steps 2.00 Z-Step Size 0 System Optimized none 😯 Z-Compensation : OFF Z Flow : 1 500 Travel Range [µm] Fast Live O Autofocus Stop

Fast Live O

10/15/2022 4:27:24 PM | Information: Stabilizing Laser Power Output.

. Choose the 'xyz' imaging mode.

2. <u>Open the Z-stack panel</u>

- Step 1. Press the 'Live' Button (1) and go to the top position of your stack using Z-focus. Click the 'Begin' button (2), which will then turn red, meaning that the 'Begin' position is locked in the top position of your stack.
- Step 2. Press the 'Live' Button (1) again and go to the bottom position of your stack using Z-focus. Click the 'End' button (3), which will then turn red, meaning that the 'End' position is locked in the bottom position of your stack.
- Now, your stack is defined.
- Define the Z-step size (4), e.g. 0.5 μm for taking the image sections every half micron.
- Click the Number of Steps option and press the Enter Key. You will be able to see the number of sections to be acquired in your stack.
- 3. To change the stack direction click the desired button under the 'Stack Direction' (5) option.
- 4. Click the start button at the bottom of the Acquire panel to image the specimen.

Take a Time-lapse



Number of Frames

Shutdown the system



- 1. Save the data and transfer the data to KCCI/data server
- 2. Stop the record time at the KCCI website.
- 3. Turn OFF the laser from the software.
- 4. Close the software and shutdown the system.
- 5. Wait.
- 6. First turn OFF the Laser button and then Power button.

*Note: Remember, do not touch the key.